



Comfort  
Release®

Global Biomedical  
Technology

# Sacral Foam Dressings for Frail or Sensitive Skin

NOW  
AVAILABLE

## 10 REASONS TO CHOOSE COMFORT RELEASE® SACRAL FOAM DRESSINGS

1. **Savings / Cost** less than silicone dressings
2. **Painless and trauma-free release** when removed as directed\*
3. **Lift & Look innovation** to check on a wound without changing the dressing\*\*
4. **Excellent moisture wicking** for moderate to heavy drainage
5. **Retains moisture** for an optimal healing environment
6. **Protective moisture-proof outer layer** to help keep out bacteria & prevent strike through
7. **Strong acrylic adhesive** for long-lasting bond to skin
8. **Low profile border** resists roll-up for longer wear time
9. **Can redistribute shear & friction** on skin
10. **Conforms / fits** gluteal cleft



Not made with natural rubber latex

*\*To remove: Simply swipe the top of the adhesive border with rubbing alcohol in prep pads or wetted cotton balls until the dressing releases from the skin easily. Rubbing alcohol not included, needs to be purchased separately.*

*\*\* Allow the alcohol to evaporate, lay the dressing back down on the skin, and the dressing will re-adhere.*

[www.comfortrelease.com](http://www.comfortrelease.com) // [info@comfortrelease.com](mailto:info@comfortrelease.com) // 833.726.9269