



Comfort Release®

Global Biomedical Technology

Next-Generation Wound Care Products

Product Catalog - Q1/2024 Update

Specially formulated for *sensitive or frail skin*.

Easy to apply, water resistant, breathable, long lasting, not made with natural rubber latex, with painless release from the skin.*

***How to Remove:** Wet the adhesive border, or surface of the tape with rubbing alcohol (alcohol prep pads or isopropyl alcohol soaked cotton balls) until product lifts easily from the skin.

Comfort Release® Adhesive Transparent Dressings

The center transparent windows are made of thin polyurethane material with a low level of adhesion, are waterproof and highly breathable for better wound healing and skin health. The border areas are made of a nonwoven film with a strong acrylic adhesive with our triggered easy release.*



- GB110** 2.375" x 2.75" Transparent Dressing (1.625" x 2" Window) | HCPCS: A6257
A Box of 10 Dressings & 20 alcohol prep pads (20 boxes per case)
- GB111** **Expected 2024:** 2" x 4" Transparent Dressing (1.25" x 3.25" Window) | HCPCS: A6257
A Box of 10 Dressings & 20 alcohol prep pads (20 boxes per case)
- GB112** 4" x 4.75" Transparent Dressing (3" x 3.75" Window) | HCPCS: A6257
A Box of 10 Dressings & 20 alcohol prep pads (20 boxes per case)
- GB112W** **Expected 2024:** 5" X 5.25" Transparent Dressing (4" x 4.25" Window) | HCPCS: A6258
A Box of 10 Dressings & 20 alcohol prep pads (20 boxes per case)
- GB113** 4.5" x 8.5" Transparent Dressing (3.5" x 7.5" Window) | HCPCS: A6258
A Box of 5 Dressings (20 boxes per case) - no alcohol prep pads included

Comfort Release® Border Adhesive Foam Dressings

Hydrophilic foam pad is made of soft foam, which is comfortable, conformable, highly absorbent and capable of expanding up to 50% when fully saturated. The breathable border has our patented easy removal adhesive.*



- GB115** 4.25" x 4.25" Border Foam Dressing (2.375" x 2.375" Foam Pad) | HCPCS: A6212
A Box of 10 Dressings & 20 alcohol prep pads (20 boxes per case)
- GB116** 8" x 8" Border Sacral Foam Dressing (6.4" x 5.4" Foam Pad) | HCPCS: A6213
A Box of 5 Dressings (15 boxes per case) - no alcohol prep pads included

Comfort Release® Adhesive Bandages

Comfort Release® Adhesive Bandages are made with a nonwoven film and a strong acrylic adhesive with our patented easy removal adhesive.*



- GB101** 1" x 3" Bandage | HCPCS: A6413
A Box w/ 20 Bandages & 20 alcohol prep pads (20 boxes per case)
- GB101-02** 1" x 3" Bandage | HCPCS: A6413
A Box w/ 100 Bandages (12 boxes per case) - no alcohol prep pads included
- GB102** 2" x 4" Bandage | HCPCS: A6413
A Box w/ 10 Bandages & 20 alcohol prep pads (20 boxes per case)
- GB102-02** 2" x 4" Bandage | HCPCS: A6413
A Box w/ 100 Bandages (10 boxes per case) - no alcohol prep pads included
- GB011** **Expected January 2024:** 1" Circular Bandage | HCPCS: A6413
A Box w/100 Bandages (48 boxes per case) - no alcohol prep pads included
- GB012** **Expected January 2024:** 0.875" x 1.25" Oval Bandages | HCPCS: A6413
A Box w/ 100 Bandages (48 boxes per case) - no alcohol prep pads included

Comfort Release® Adhesive Tapes

Comfort Release® Adhesive Tapes are made with a transparent, bidirectional tear, polypropylene film. Made with a strong acrylic adhesive and our patented easy removal adhesive.*



- GB103** 1" x 4 Yard Transparent Tape | HCPCS: A4450
A Box of 12 Rolls (10 boxes per case) no alcohol prep pads included
- GB103-RET** 1" x 4 Y Transparent Tape | HCPCS: A4450
One Role of Tape (24 rolls per case) - no alcohol perp pads included
- GB104** 2" x 4 Yard Transparent Tape | HCPCS: A4450
A Box of 6 Rolls (10 boxes per case) - no alcohol prep pads included

Comfort Release® Supplies

Additional wound care supplies.

- 861** Dukal Alcohol Prep Pads Box of 100 Large Alcohol Prep Pads (10 boxes per case)
- 856** Dukal Sting Relief Pads Box of 200 Sting Relief Pads (20 boxes per case)
- 887** Alcohol Swabsticks - Bag of 15 Antiseptic 4" Swabsticks, individually packaged

Skintact® ECG Electrode Pads

Made with a nonwoven polyurethane film border, water resistant, breathable, long lasting with painless, trauma-free removal. Lifts off without pulling hair. Made for cardiac monitoring and diagnostic procedures. Anyone with sensitive skin.

- GB135** **Expected January 2024:** 2" Diameter, round, dry gel with Comfort Release® surrounding border - A Pouch of 15 Pads



Comfort Release® Surgical Drape

Strong medical grade acrylic ensures a strong and effective seal while Comfort Release® patented adhesive allows the drape to be removed painlessly and trauma-free. The drape acts as a barrier to external contaminants and provides waterproof protection.

- GB141-01** **Expected 2024:** 10.5" x 12" Drape - A Pouch with 1 Drape



How to Apply and Remove Comfort Release® Products

To Apply: Wash your hands. Clean the wound and the skin around the wound. Allow the skin to dry. Apply the product to the skin, pressing firmly to ensure adherence.

To Remove: Use alcohol prep pads or cotton balls soaked with rubbing/isopropyl alcohol and swipe across the surface of a bandage or tape or across the border of a dressing. Wait until the product is saturated, then lift easily from the skin.



To order Comfort Release® Wound Care Products, go to: [comfortrelease.com/where-to-buy](https://www.comfortrelease.com/where-to-buy) or call **470.636.4410**. For more info, email us at info@comfortrelease.com