



Comfort Release®

Global Biomedical Technology

Next-Generation Wound Care Comfort Release® Tapes

TAPE THAT CAN BE APPLIED TO THE SKIN WITHOUT FEAR OF INJURY

TO REMOVE PAINLESSLY: Wet the top surface of the tape with rubbing alcohol (prep pads) until the tape lifts easily from the skin.



Avoid skin stripping, skin tears, rashes, erythema, contact dermatitis and maceration.



1" & 2" Wide Perforated Tape

- Tears easily without scissors
- Water resistant
- Breathable
- Long-lasting
- Hypoallergenic
- Not made with natural rubber latex

